

B/W Digital Trolley Ultrasound Workstation System



We supply **Clinical Equipment**(Color ultrasound, B/W ultrasound, Human ultrasound, Veterinary ultrasound, Portable ultrasound, Trolley ultrasound, etc...),we devoted ourselves to send you clients high quality Clinical Equipment with long terms services.

We are expecting to become your long-terms business partner in China.

Introduction of Digital Trolley Ultrasound Scanner with Workstation YJ-U80T:

Color doppler YJ-U80T is a Digital Trolley B/W Ultrasound Diagnostic System machine, equipped with, 15" LCD monitor. The device is applicable for hospitals and clinics for examining liver, GB, spleen, kidney, pancreas, heart, bladder, uterus, etc.

Note:

1,Two years quality guarantee and 24 hours phone service for all the products supplied.

2, Training Service Supplied!

Specifications of Full Digital Ultrasound Scanner YJ-U80T:

Full Digital Ultrasound Scanner YJ-U80T	
Ultrasound Scanner	
Detecting depth	≥ 220mm
Monitor	15" LCD monitor
Image Modes	B, B+B, B+M, M, 4B
Gray Scale	256
Image amplification	×1.0; ×1.2; ×1.5; ×2.0
Local zoom	2 times
Electronic focusing	4 focuses combination randomly
Measurement	distance; circumference; area; heart rate; gestational week; fetal weight etc
Annotation	time; date; ID; age; sex; detecting depth; probe type; focus; frame correlation etc
Cine-loop	341
Image load/storage	32
Body marks	40 types
Obstetric table	BPD; CRL; GS; FL; AC; HC; FW; EDD
Language	English-Chinese
Power requirements	100V~130V/220V~240V, 50/60Hz
Net Weight	40Kgs(including probe)
Dimension (mm)	910*760*900; 40*35*37 (CRT monitor carton)
Imaging workstation specification	
Main board	Gigabyte G31M
CPU	INTEL Celeron 430 1.8MHz
EMS memory	Kingston800 1G
Hard disk	Western Digital 160G
Monitor	19" LCD
Electromechanica l	phillips+ GREATWALL
Keyboard	Lisheng
Mouse	double swallow optional mouse
Standard configuration:	1) Ultrasound Scanner with 17" CRT Monitor: 1pc 2) Image Workstation: 1pc

	3) 3.5MHz Convex Probe: 1pc 4) 6.5MHz Trans-vaginal Probe: 1pc
Optional	7.5 MHz Linear Probe Video printer and Ink-jet printer
HS Code	901812910

Features of Full Digital Ultrasound Scanner YJ-U80T:

- Display image in real time, frozen, amplified, black- white reversed, up/down conversed and shift depth in real time.
- Video output: PAL and VGA; large screen monitor, videocassette recorder, video image printer.
- Soft keyboard and track ball to make operation faster, convenient and flexible.
- Trolley structured with plastic cover. Switch power is used without working frequency transformer.
- FPGA and SMT in the device make the device high integrated, small volume and light weight.
- The device passed clinical verifications of safety and diagnosis validity in its manufacturing

FAQ:

Q:What is time of delivery?

A:Since our products are produced according to each specific order, the delivery period will be 15days after confirmation of order. If it is small or sample order, please do check with us about whether we have ready goods in stock to make delivery as soon as possible.

Q:What about quality assurance?

We offer one year warranty period beginning from the time of shipment for medical equipment, Per minor problems, we can provide free spare parts for replacements, Per serious problems, we can make replacement of free charge. But if the machine is damaged due to rough handling during the transportation, we are free of any responsibility, but we will still provide help to the buyer issue claims against the shipping company or insurance company, thus all the spare parts or expenses should be borne by buyer.

Q:What is payment terms?

A:T/T in advance, Western Union

Q:Which shipping company usually using?

A:MAERSK-LINE, PIL, COSCO, SAFMARINE, MARUBA, APL, CMA, CSAV, K-LINE, MSC, NYK, AND DHL, EMS, TNT, UPS FOR COURIER, buyer's option will be our priority, we will give suggestion according to the schedule and services.

Delivery:



Related products:

(YJ-U80TPLUS)

